LAICA®

LAICA S.p.A.
Viale del Lavoro, 10 - Fraz. Ponte
36021 Barbarano Vicentino (VI) - Italy
Tel. +39 0444.795314 - 795321
Fax. +39 0444.795324
Made in China

www.laica.com

LAICA®

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A new electronic personal scale that allows to evaluate your body composition and to automatically calculate the body fat and water percentage in function of your height, age and sex.

The exact knowledge of your body composition (calculation of lean body mass, fat body mass and total body water contents), allows a correct estimate of the excess fat and at the same time determines the advisable weight as a function of your nutritional condition.

Furthermore, the knowledge of your body composition allows to combine a correct nutrition and physical activity as essential conditions to improve the quality of life, promote health and prevent pathologies. It has a memory function with a 10-person capacity.

Capacity 180 kg / 396 lb / 28 st 4 lb, divisione 100 g / 0.1 lb / 0.1 lb (< 10 st) 1 lb (> 10 st)

Warning!

- **The body fat and water percentage is to be considered only as approximate information.**
- **For any further information, please consult your physician or dietician.**
- **Self-measuring means control, it does not mean diagnosis or treatment.**
- **Any unusual levels should always be discussed with your doctor.**
- **Under no circumstances must you change the dose of medicines prescribed by your Doctor.**

IT IS IMPORTANT TO CAREFULLY READ THE INSTRUCTIONS AND WARNINGS IN THIS BOOKLET BEFORE USE AND TO KEEP THE BOOKLET IN A SAFE PLACE.

BODY COMPOSITION

It is important to know the body composition of an individual in order to be able to estimate his/her health condition.

Actually, it is not enough to know one's weight but it is necessary to define its quality, which means dividing the human body in its two major component masses:

- lean body mass: bones, muscles, viscera, water, organs, blood;
- fat body mass: adipose cells.

Under the same conditions of weight, these two masses have different volumes.

By means of electrodes placed on the platform, the scale sends a weak low intensity electric current that passes through the body mass.

The electric current reacts in a different way according to the masses it passes through: the lean body mass opposes a weak resistance, while the fat body mass has a greater resistance.

This body resistance, called bioelectrical impedance, varies also in function of sex, age and height. For a correct control of one's ideal weight, it is important, besides the weight, to compare the body fat percentage calculated by the scale, with the parameters contained in the table and the body water percentage, with the information listed hereinafter.

If the percentages are outside the normal values, we suggest you to consult your doctor.

In order to better understand the value of your body fat percentage, we suggest you write down this value for a given period of time. In order to obtain a valid comparison, you should weigh yourself

always at the same time of the day and under the same conditions.

The best moment to weigh yourself is in the morning, before breakfast, without clothes. You have to wait 15 minutes after getting up, so that the water get distributed itself throughout the body.

Body fat mass percentage

The percentage of fat listed in the following table, represents the total fat, calculated by summing up the essential fat (or primary fat) and the spare fat (or storage fat).

The essential fat is necessary since it is needed for the metabolism. The percentage of this fat is different between the two sexes: about 4% for men and about 12% for women (percentage calculated on the total body weight).

The spare fat, which is theoretically superfluous, is necessary in a minimum quantity for a healthy condition, especially for those who carry out physical activities.

Even in this case, the fat percentage is different between the two sexes: 12% for men, 15% for women (percentage calculated on the total body weight).

Therefore, the normal percentage of fat mass on the body weight is approximately 16% for men and 27% for women. Such percentages differ with the person's age as shown in the table here below.

BODYFAT TABLE					
	Age	Excellent	Good	Fit	Overweight
MAN	15-39	<11%	11.1-23%	23.1-31%	>31.1%
	40-59	<12%	12.1-24%	24.1-32%	>32.1%
	60-79	<13%	13.1-25%	25.1-33%	>33.1%
WOMAN	15-39	<20%	20.1-32%	32.1-38%	>38.1%
	40-59	<21%	21.1-33%	33.1-39%	>39.1%
	60-79	<22%	22.1-34%	34.1-40%	>40.1%

Water percentage

Water is the main component of lean body mass; the water percentage over the fat body mass is 73%. The lean body mass is obtained by subtracting from the total weight the fat body mass (for instance for a 25-year old man, with a fat body mass of 17%, lean body mass is equivalent to 83%).

About 60% of the body weight is made up by water: it therefore represents the fundamental component of every living being and it is necessary to keep us alive.

Women (having normally more adipose tissue than men), have a water percentage of about 55-58%, while men have about 60-62% of it.

The water percentage is greater in childhood (at birth it varies from 70 to 77%), and it decreases with age (in old people the percentage varies from 45% to 55%) and with the increase of fat deposits.

Obese people have a lower body water percentage than lean people.

Refer to this data to compare the obtained percentage.

Instructions and guarantee

Water accomplishes many biological functions:

- it promotes the digestive processes;
- it is a thinner for substances that are important for the organism;
- it removes metabolic wastes;
- it carries oxygen, nutrients, enzymes, hormones and glucose to the cells;
- it is a source of mineral salts;
- it regulates the body temperature.

The human organism excretes every day a given quantity of water, which must be regularly replaced; a water loss of approximately 10% of the body weight, may put at risk the health condition. Body weight variations, may occur throughout short (hours), medium (days) or long (weeks) periods.

Small weight variations that take place in a few hours, are completely normal and are due to water retention, since the extra-cellular water (which includes interstitial fluid, plasma, lymph and transcellular fluid) is the only compartment subject to such rapid variations.

When a weight reduction occurs in a few days, it is very likely that it affected the total water amount (distributed by approximately 60% inside the cells - intracellular water - and by approximately 40% outside the cells - extracellular water-) and therefore, it also affected the lean body mass. To conclude, the fat body mass is subject to variations which can be recorded in weeks.

Therefore, in order to evaluate the progression of one's body composition it is necessary to keep in mind that:

- **a weight reduction by itself, accompanied by an unchanged or even increased body fat percentage, shows that the body has simply lost water (for instance after a training session, a sauna or a diet that is limited to rapid weight loss).**
- **a weight increase accompanied by a decreasing or unchanged body fat percentage, shows the development of precious muscular mass.**
- **when weight and body fat decrease at the same time, it means that the diet is effective and that we are also losing volume.**
- **when there is a high percentage of fat or a reduced muscular mass, the body water percentage is always too low.**

Warning!

When following a diet, it is necessary to expect a reduction in the fat body mass and not in the muscular mass.

Examples of body composition

Person 1

Entered data: man, 185 cm tall, 39 years old.

Data calculated by the scale: Weight 105 kg , Bodyfat 30.3%, Bodywater 51%

Analysis of body composition: An overweight body composition results from the table. Let's study in any case the analysis in deep.

The percentage of lean mass is 69.7% (100%-30.3%). Body fat consists of:

- essential fat => 105 kg x 4%= 4.2 kg
- reserve fat => 105 kg x 12%= 12.6 kg

Total body fat (essential fat + reserve fat) is 4.2+12.6 kg = 16.8 kg, which compared to the total body weight (105 kg) is equal to 16% (ideal % of fat mass for person 1).

The difference between 30.3% (actual % of fat mass) and 16% (ideal % of fat mass) is 14.3%, which applied to the total weight (105 kg) is equal to 15 kg.

15 kg represent the kg in defect of person 1.

The water percentage of 51% is low compared to the normal percentage that must have the men but it is due to the high percentage of fat mass and then a reduced muscular mass.

TIPS FOR USE

- **Place the scale on a flat, hard surface.**

Soft, uneven flooring (e.g. rugs, carpets, linoleum) can cause the scale to give an uneven reading of your weight.

- Step onto the scale with your feet parallel and weight equally distributed (see illustrations). Stand still while the scale measures your weight.
- It is important to weight yourself with dry, bare feet, correctly positioned on the electrodes (see illustrations). Moreover, for the scale to work properly, avoid your legs touching.
- **The scale does not show any danger, since it works with 4 x 1,5V alkaline batteries. However, the following people should avoid using it:**
 - those fitted with a pace-maker or other forms of electronic, medical equipment
 - those with a temperature
 - children under the age of 10
 - pregnant woman
 - those in dialysis, with edemas on their legs, people with dimorphism or suffering from osteoporosis.
 - people on cardiovascular medication.

Always weigh yourself on the same scale every day at the same time, undressed and before breakfast.

To get the best results from your scale, weight yourself twice, and if the two weights are different from each other, your weight is between the two readings.

Wait for about 15 minutes after getting up, so that the water can distribute through around your body.

- The surface of the scale can be slippery, if wet.
- Warning!
- The measurement function of body fat and water percentage, may be utilized only by people that belong to an age group that goes from 10 to 80 and that are between 100 and 240 cm tall.
- The precision of this fat and water percentage measurement can be altered in the following cases:
 - a) a high alcohol level
 - b) a high level of caffeine or drugs
 - c) after an intense sports activity
 - d) during an illness
 - e) during pregnancy
 - f) after heavy meals
 - g) in presence of dehydration problems
 - h) your personal data (height, age, sex) have not been input correctly

- i) if you have wet or dirty feet
j) people with swollen legs

Warning!

The body fat and water percentage is to be considered only as approximate information.

For any further information, please consult your doctor or dietician.

OPERATING INSTRUCTIONS

a. How to insert batteries

This scale operates with 4 replaceable 1.5 V alkaline batteries. When the batteries are dead the display shows "LO" or the battery symbol. Insert new batteries with the positive pole upwards. Remove the run down batteries by lifting the cover of the battery compartment located on the bottom of the device.

Extract the batteries and dispose them as special waste at the specific collection points forecast for recycling.

For more information about the disposal of run down batteries, please contact either the store where the device containing the batteries was purchased, or the Municipality, or the waste disposal local service.

Caution: For greater safety, children younger than 12 years should not remove the batteries.

b. Unit of measure key

With the scale switched off, select the measuring unit (kg, lb or st) via the key on the bottom of the scale.

c. Use as simple scale

- 1) Press "ON" button to turn the scale on.
- 2) Wait a few seconds and when "0.0" appears on the display screen, the scale is ready for use.
- 3) Step onto the scale and wait without moving; after a few seconds, the display screen will show your weight. Do not lean on anything while being weighed as this could alter the reading.
- 4) Step off the scale, the display screen will continue to show your weight for a few seconds then it will turn itself off automatically.

Each time the scale is moved and then returned to its resting surface, the first weight must not be considered due to reorganization of the electronic board with the loading cells. Consider the second and successive weights as in certain cases the first one could be altered for the reason mentioned above.

d. How to input personal data

This scale may memorize the data of 10 different persons.

- 1) Press the "SET" button.
- 2) "P1" will flash on the display screen, indicating the "01" person: if you wish to change the memory location (memories from 0 to 9), press the "▲" or "▼" symbols. Once you have selected the memory location, confirm by pressing "SET".

- 3) Continue by pressing the "▲" or "▼" symbols to define your sex. Confirm by pressing "SET".
- 4) Use the "▲" or "▼" symbols again to insert age (from 10 to 80) and then confirm by pressing "SET".
- 5) Insert your height (from 100 cm to 240 cm) using the "▲" or "▼" symbols to increase or decrease the value. The display screen will give the value in cm or ft. Press SET to confirm.
- 6) Now the scale has stored all your data; the display shows them to allow you checking their correct entering.
At this point, the display shows "0.0". Get on the scale paying attention to position the naked feet correctly on the two electrodes. Then, the display shows the following data sequentially:
- your weight
- the message FAT (bodyfat) and the percentage of fat present in the body
- the message TBW (bodywater) and the percentage of water present in your body
Compare the obtained fat percentage (FAT) obtained with the parameters given in the "Bodyfat table": if the values are not included in the normal range, refer to your doctor.
Estimate the water percentage (TBW) as in paragraph "% of water".
- 7) The scale switches off automatically after about 10 seconds.

To modify the entered data, call up the involved memory area and repeat the whole proceeding. To enter the data of another person, select a free memory area and repeat the whole proceeding.

Caution: while being in the "How to input personal data" function, the scale shall switch off if not any key is pressed within 10 seconds.

e. Calculation of weight %, of body fat % and of body water %

- 1) Activate the scale by keeping pressed button ON until the desired memory zone is shown. Previously entered data will be displayed by the display, followed by "0.0kg".
- 2) Get on the scale and await the display of your data (see item 6 of the paragraph "introduction of personal data").

Caution!

If the display does not show "0.0", it will not be possible to read the weight, as well as to calculate the percentage of fat and water. If the display shows "Err", it means that an error occurred during the measurement. Repeat the whole procedure again. In case of overweight, the display shows "EEEE". In this case, remove immediately the weight from the balance in order not to damage its mechanism.

TECHNICAL FEATURES

Capacity: 180 kg / 396 lb / 28 st 4 lb

Division: 100 g / 0.1 lb / 0.1 lb < 10 st / 1 lb > 10 st

Measurement unit: kg/lb/st

Memories for 10 persons

Four 1,5V alkaline batteries type AAA

Accuracy: +/- 1% + 1 division of applied weight

LCD display

Storage environmental requirements: +10°C + 50°C; RH ≤85%

Operating environmental requirements: +10°C + 40°C; RH ≤85%

Instructions and guarantee

MAINTENANCE

- Handle your scale with care as it is a precision instrument. Therefore, avoid to expose it to extreme temperatures, damp, knocks, dust and direct sunlight. Do not drop the scale or knock on it.
- Keep the scale away from sources of heat.
- **Do not jump on the scale.**
- **Clean the scale with just a soft dry cloth.**
- **To clean the scale never use chemical or abrasive products.**
- Take care that no liquids penetrate into the scale, as it is not waterproof. Do not leave or use the scale in very damp surroundings, humidity above 85% or contact with water and other liquid parts could cause corrosion and jeopardize the precision of the scale.

Important! Read the instructions carefully

This product is meant for domestic use.

The conformity, certified by the CE-marking on the device, refers to the Directive 30/2014/Eu on electromagnetic compatibility.

This product is an electronic device that was tested to assure, according to the current technical knowledge, that it does not interfere with other nearby located equipment (electromagnetic compatibility) and that it is safe if used according to the advices given in the instructions for use. In case of operation anomalies, do not keep using the device and, if necessary, refer to the manufacturer directly.

DISPOSAL PROCEDURE (Dir. 2012/19/Eu-WEEE)

The symbol on the bottom of the device indicates the separated collection of electric and electronic equipment.

At the end of life of the device, do not dispose it as mixed solid municipal waste, but dispose it referring to a specific collection centre located in your area or returning it to the distributor, when buying a new device of the same type to be used with the same functions.

If the appliance to be disposed of is less than 25 cm, it can be returned to a retail location that is over 400 m² without having to purchase a new, similar device.

This procedure of separated collection of electric and electronic devices is carried out forecasting a European environmental policy aiming at safeguarding, protecting and improving environment quality, as well as avoiding potential effects on human health due to the presence of hazardous substances in such equipment or to an improper use of the same or of parts of the same.

Caution! The wrong disposal of electric and electronic equipment may involve sanctions.

WARRANTY

This equipment is covered by a 2-year warranty starting from the purchase date, which has to be proven by means of the **stamp or signature of the dealer and by the fiscal payment receipt that shall be kept enclosed to this document.** This warranty period complies with the laws in force and applies only in case the consumer is an individual.

Laica products are designed for domestic use only; therefore, their use for commercial purposes is forbidden.

The warranty covers only production defects and does not apply, if the damage is caused by accidental events, misuse, negligence or improper use of the product. Use the supplied accessories only; the use of other accessories invalidate the warranty cover.

Do not open the equipment for any reason whatsoever; the opening or the tampering of the equipment, definitely invalidate the warranty. Warranty does not apply to parts subject to tear and wear along the use, as well as to batteries, when they are included.

Warranty decays after 2 years from the purchase date. In this case, technical service interventions will be carried out against payment. Information on technical service interventions, under warranty terms or against payment, may be requested by writing at info@laica.com.

Repair or replacement of products covered by the warranty terms are free of charge. In case of faults, refer to the dealer; DO NOT send anything directly to LAICA.

All interventions covered by warranty (including product replacement or the replacement of any of its parts) shall not be meant to extend the duration of the original warranty period of the replaced product.

The manufacturer cannot be made liable for any damages that may be caused, directly or indirectly, to people, things and pets resulting from the inobservance of the provisions contained in the suitable user manual, especially with reference to the instructions related to installation, use and maintenance of the equipment.

Being constantly involved in improving its own products, Laica reserves the right to modify its product, totally or partially, with reference to any production need, without that this involves any responsibility whatsoever by Laica or its dealers.

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