IT Italiano

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ELECTRONIC PERSONAL SCALE

Instructions and warranty

Your new bathroom scale is extremely reliable and precise, it does not require any sort of maintenance. To improve the reading, weigh yourself twice, and if the weights are different, your weight is the middle of the two readings. Always weigh yourself on the same scale, at the same time of day, preferably undressed and before breakfast. Attention! The surface of the scale can be slippery, if wet.

IT IS IMPORTANT TO READ THE INSTRUCTIONS AND WARNINGS IN THIS HANDBOOK BEFORE USE AND KEEP IT IN A SAFE PLACE.

THE BODY WEIGHT

Checking the weight means performing an **energetic balance**: such balance results from the relation between the quantity of ingested nutrients, i.e. the ingested calories, and the actually consumed calories. The daily necessary quantity of energy depends on the metabolism, on the life-style and on the type of activities performed by each single person.

Controlling one's own weight is essential in order to protect health effectively: eating well means being healthy. It has been proved that most of diseases, many of which are extremely dangerous, are associated to unbalanced eating habits not only in terms of quality but also in terms of quantity.

- Among the most known, the following ones should be mentioned:
 - ischemic heart diseases;
 - cerebrovascular diseases;
 - diabetes mellitus;
 - intestine tumors;

All these diseases are present above all among obese persons; therefore, controlling the weight together with controlling the quality of the nutrients becomes one of the essential factors to prevent diseases.

Healthy eating means selecting the foodstuffs that supply the body with all the necessary nutriment to grow and keep the physiological functions and the biochemical processes balanced.

Hints for a healthy feeding and to control weight:

- To follow a balanced diet it is necessary not to exceed in the consumption of the same types of food neglecting the essential elements contained in other types of food, as well as not eating more than needed.
 - A healthy diet includes suitable quantities of proteins, carbohydrates, fats, fibers, vitamins and minerals, as well as, obviously, water.
- Have a meal at regular intervals. It is a bad habit. It is better to have many small meals instead of just an abundant one.
 - Three meals a day and a sneak in the middle of the morning or of the afternoon represent the basis for a correct feeding.
- It is suggested to sit down and relax when eating. Select a comfortable environment can be relevant.
- Do not eat when you are focused on doing other activities, i.e. when you are watching TV or working, in order to avoid eating too much without even noticing it.
- To support the digestion process it is recommended to eat slowly and chew well, enjoying the color, taste and consistency of the food.
- Combining a correct feeding to physical activity is an essential condition to allow increasing the metabolism, consuming more calories and remaining in good shape.
- It is important to remember that with age human body seems to need less energy due to still not so clear reasons (maybe due to the fact that one leads a less active life); however, many aged people do not reduce their consumption of food and tend to grow fat.
 - Getting used to follow a balanced diet is surely of help in increasing the quality of life, promoting health and preventing diseases, as well as in increasing vitality and remaining in good shape!

Caution!!!

- Never follow diets autonomously, always refer to a doctor or dietician.**
- Self-measurement means control, not diagnosis or treatment.**
- Unusual values have always to be discussed with one's own physician.**
- Under no circumstances must you change the dose of medicines prescribed by your Doctor.**

I NSTRUCTIONS FOR USE

- Select the weighing unit of measure (kg, lb or st) with the key located on the bottom of the scale.
- It is possible to select the unit of measure if the scale is turned on.
- Place the scale on a flat, hard surface. If the floor is loose or uneven (e.g. mats, carpets, linoleum) the exactness of your weight is not guaranteed.**
- Step onto the scale with your feet parallel and your weight equally distributed (see illustrations).**

Warning: the scale switches on once step on. Stand still while the scale measures your weight. Do not lean on anything while being weighed as this could alter the reading.
- Step off the scale, the display screen will continue to show your weight for a few seconds then the scale will turn itself off automatically.

WARNING:

- When the display screen shows "LO" or the symbol of the battery, it's necessary to replace the batteries bearing in mind the polarity indicated (+) / (-).
 - In case of overweight, the display screen shows "Er", "--" or "OL", remove the weight from the scale to avoid any damage of the mechanism.
- BATTERIES**
- This scale operates with one 3V replaceable lithium battery, type CR2032. Remove the protective film from the batteries before use. Remove the battery if the scale is not to be used for a prolonged period of time.
- extract the screws from the cover of the battery compartment.
 - open the cover and remove the battery
 - insert the new battery and place the positive pole upwards.
 - place the cover back on the battery compartment.
 - screw the screws in the cover.

Extract the battery and dispose it as special waste at the specific collection points forecast for recycling. For more information about the disposal of turn down batteries, please contact either the store where the device containing the batteries was purchased, or the Municipality, or the waste disposal local service.

Caution: For greater safety, children younger than 12 years should not remove the batteries.

MAINTENANCE

- Handle your scale with care as it is a precision instrument.
- The scale should avoid to expose to extreme temperatures, damp, knocks, dust and direct sunlight.
- Do not drop the scale or knock on it.
- Keep the scale away from sources of heat.
- Do not jump on the scale.**
- Clean the platform with a damp cloth.**
- Do not clean the bottom of the scale with wet or soft dry cloth.**
- Do not use the scale with any chemical or abrasive products.**
- Take care that no liquids penetrate into the scale, as it is not waterproof.
- Do not leave or use the scale in very damp surroundings, humidity above 85% or contact with water and other liquid parts could cause corrosion and jeopardize the precision of the scale.

TECHNICAL DATAS

One 3V CR2032 lithium battery

LCD display

Maximum tolerance: +/-1% + 1 division

Working conditions: +10°C +40°C; RH 85%

Storage conditions: +10°C +50°C; RH 85%

i Important! Read the instructions carefully

GB English

ELECTRONIC PERSONAL SCALE

Instructions and warranty

This scale is meant for domestic use.

The conformity, certified by the CE-marking on the device, refers to the EEC Directive 2004/108 on electromagnetic compatibility.

This product is an electronic device that was tested to assure, according to the current technical knowledge, that it does not interfere with other nearby located equipment (electromagnetic compatibility) and that it is safe if used according to the advices given in the instructions for use.

In case of operation anomalies, do not keep using the device and, if necessary, refer to the manufacturer directly.

DISPOSAL PROCEDURE (Dir. 2012/19/Eu-WEEE)

The symbol on the bottom of the device indicates the separated collection of electric and electronic equipment.

At the end of life of the device, do not dispose it as mixed solid municipal waste, but dispose it referring to a specific collection centre located in your area or returning it to the distributor, when buying a new device of the same type to be used with the same functions.

If the appliance to be disposed of is less than 25 cm

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