LAICA[®]

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STANDARDS

Il prodotto risponde ai seguenti standards: IEC 60601-1 (Medical electrical equipment - Part 1:General requirements for basic safety and essential performance), IEC 60601-1-2 (Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests), IEC 80601-2-30 (Medical electrical equipment – Part 2-30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers).

Prodotto da: Shenzhen Pango Medical Electronics Co., Ltd. No. 25 1st Industry Zone, Fenghuang Rd, Xikeng village, Henggang Town, Longgang District, Shenzhen, Guangdong China

Lotus NL B.V., Koningin Julianaplein 10, 1e Verd, 2595AA, The Hague, Netherlands

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ARM BLOOD PRESSURE MONITOR TYPE PG-800B10
INSTRUCTIONS AND WARRANTY

Dear customer, Laica would like to thank you for choosing one of our products, designed according to strict performance and quality criteria in order to ensure total satisfaction.

IMPORTANT
READ CAREFULLY BEFORE USE
KEEP IN A SAFE PLACE FOR FUTURE REFERENCE

The instruction manual must be considered an integral part of the product and must be kept for its entire life. If the product is transferred to another owner, its documentation must also be transferred in its entirety.
To ensure the safe and correct use of the product, the user must carefully read the instructions and warnings contained in the manual insofar as they provide important information regarding safety, user and maintenance instructions.
Should the instruction manual get mislaid or you require additional information or clarifications, please fill in the relevant form on the website: <https://www.laica.it/> in the FAQ and Support section.

This completely automatic and non-invasive device measures and checks blood pressure (systolic and diastolic), the heart rate and checks for arrhythmias.

CONTENTS	
KEY TO SYMBOLS	pag. 12
SAFETY WARNINGS	pag. 12
CLASSIFICATION OF BLOOD PRESSURE READINGS	pag. 13
PRODUCT DESCRIPTION	pag. 14
INSTRUCTIONS FOR USE	pag. 14
MAINTENANCE	pag. 16
TROUBLESHOOTING	pag. 16
DISPOSAL PROCEDURE	pag. 16
WARRANTY	pag. 17
STANDARDS	pag. 17
ELECTROMAGNETIC COMPATIBILITY	pag. 18

Warning

Prohibition

Caution! Read the user instructions carefully

Symbol for "Type BF applied parts" (the cuff is type BF applied part)

Keep dry!

C

CE

0413

EC

REP

LOT

In compliance with European legislation on medical devices

European Representation

Batch production number

Manufacturer

SAFETY WARNINGS

- IP20:** Degree of protection of electrical equipment covers, where the first figure indicates the degree of protection against penetration by solid foreign bodies (from 0 to 6) and the second figure the degree of protection against penetration by liquids (from 0 to 8).
- Prior to use, ensure the product is intact and that there is no visible damage. If in doubt, refrain from using the product and contact the retailer that sold you the product.
 - Keep the plastic bag away from children: danger of suffocation.
 - This device must be used only for its intended purpose and according to the user instructions. All other types of use shall be considered inappropriate and therefore hazardous. The manufacturer shall not be deemed liable for any damage caused by inappropriate or incorrect use.
 - People with physical, sensory or mental disabilities or unskilled people, may use and perform maintenance on the product only under adult supervision. Children must not play with the device.
 - To prevent the possibility of accidental strangulation, keep the device away from children and avoid wrapping the cuff around the neck.
 - Handle the product with care, protect it against accidental impact, extreme temperature fluctuations, humidity, dust, direct sunlight and sources of heat.
 - In the event of a failure and/or malfunction, turn off the device without tampering with it. Always contact your retailer for repairs.

CAUTION! PRIOR TO USING THIS DEVICE

- Ensure hands are dry when using the on and off key.
- ⚠ NEVER submerge the product in water or any other liquids.

- This device can measure blood pressure in adults, 18 years and older (arm circumference from approx. 22 to 42 cm). To measure a child's blood pressure, consult your physician.
- ⚠ DO NOT use this device if you are affected by severe arrhythmias.
- Measuring your own blood pressure is a way of monitoring it and must not be considered a diagnosis or treatment. Unusual readings must always be discussed with your physician. Under no circumstances should you alter the dosages of any medication prescribed by your physician.
- Consult your physician prior to using the device in the following cases:
 - patients with heart pacemakers,
 - patients with an irregular heartbeat (arrhythmia),
 - pregnant women,
 - when applying the device on an arm wound or injury,
 - when applying the device on a limb where there may be an intravascular access device or arteriovenous shunt (A-V),
 - when applying the device on a person that has had a mastectomy,
 - when using the blood pressure monitor at the same time as other medical monitoring devices already being used on the same limb,
 - when undergoing dialysis treatment,
 - when taking anticoagulants, antiplatelets or steroids.
- In the following cases, the reading may be lower or imprecise: arteriosclerosis, upper limb spasticity, low blood circulation, cardiovascular system disorders, very low blood pressure, blood flow disorders, arrhythmias and other pre-pathological conditions.
- The device must not be used for first aid and for the continuous monitoring of blood pressure.
- This device may only be used in domestic environments, do not use it outside or while travelling.
- The device may give inaccurate readings if it is used in temperatures or humidity levels outside the limits indicated in the "Technical Specifications" paragraph.
- ⚠ Do not use near strong magnetic fields and therefore keep away from radio or mobile phone systems (for more information on interferences, see the paragraph "Electromagnetic compatibility").

- Use only with the original manufacturer's cuff. Use of non-original cuffs may lead to incorrect readings.
- ⚠ Do NOT share the use of the cuff with persons with infectious diseases (risk of infection).
- ⚠ Do NOT bend or pinch the air tube during the measurements so as to avoid possible errors of cuff inflation or bruising of the arm due to the continuous pressure in the cuff.
- SAFE USE OF BATTERIES**
 - Remove the batteries if the device is not used for long periods of time and keep them in a well-ventilated and dry place, at room temperature.
 - ⚠ Do NOT recharge the batteries if they are not rechargeable.
 - ⚠ Do NOT charge rechargeable batteries using methods other than those indicated in the manual or using unsuitable equipment.
 - ⚠ NEVER expose the batteries to sources of heat or direct sunlight. Failure to observe these instructions may result in damage to and/or the explosion of the batteries.
 - ⚠ Do NOT throw the batteries into an open flame.
 - Batteries must always be removed or replaced by adults only.
 - Keep batteries out of reach of children: the ingestion of batteries can cause serious harm and death. If a battery is swallowed, seek urgent medical care.
 - The acid in batteries is corrosive. Avoid contact with the skin, eyes or garments.

CLASSIFICATION OF BLOOD PRESSURE READINGS			
BLOOD PRESSURE CLASSIFICATION	SYSTOLIC (mmHg)	DIASTOLIC (mmHg)	SEGMENT COLOR
Optimal	<120	<80	Green
Normal	120 – 129	80 – 84	Green
Normal – High	130 – 139	85 – 89	Yellow
Ipertension 1st grade - Light	140 – 159	90 – 99	Orange
Ipertension 2nd grade - Moderate	160 – 179	100 – 109	Orange
Ipertension 3rd grade - Severe	≥180	≥110	Orange

Values below 105 mmHg (systolic) and 60 mmHg (diastolic) indicate hypotension. We recommend seeking medical advice.
This device is able to detect irregular heartbeats, or arrhythmia, indicating them on the display with the symbol (♥).
Arrhythmia may be caused by frequent anxiety, certain emotional states, excessive alcohol use, genetic predisposition, age or other.
It can be a symptom of a particular physical or mental condition (temporary disturbance) or an actual heart problem.
Always seek medical advice if the monitor displays the irregular heartbeat symbol.

PRODUCT DESCRIPTION (see fig.1)

- 1) LCD DISPLAY
- 2) "SET" key
- 3) "M" key
- 4) "O/I" key
- 5) Cuff storage compartment
- 6) Air hose socket
- 7) Battery compartment
- 8) Cuff
- 9) Air hose
- 10) Air hose plug
- 11) Batteries
- 12) Case

TECHNICAL SPECIFICATIONS

- Product name: Blood pressure monitor
- Commercial name: BM2605
- Classification: Internally powered, Type BF applied part, IP20, No AP or APG, Continuous operating mode
- Method: oscillometric, automatic air inflation and measurement
- Measurement range: from 30 to 280 mmHg (cuff pressure), from 40 to 199 beats/min (heart rate)
- Precision: systolic, diastolic ±3 mmHg heart rate ±5% of reading
- Memories: 90 memories
- Arm circumference: between approx. 220 and 420 mm

- Power supply: 4 x 1.5V AA alkaline batteries (LR06)
- Battery autonomy: approx. 200 readings
- Environmental conditions for use: from +5°C to +40°C; relative humidity ≤15~93% RH
- Environmental conditions for storage: from -20°C to +55°C; relative humidity ≤93% RH
- Atmospheric pressure for use and storage: 50kPa-106kPa
- Life expectancy of product with use: 5 years

INSTRUCTIONS FOR USE

BATTERY INSERTION/REPLACEMENT

This blood pressure monitor works with 4 alkaline batteries of the 1.5V AA type.
At the first use, and when the display shows the battery symbol (🔋), proceed with the insertion and /or replacement of the batteries.
Open the battery compartment by sliding the lid in the direction of the arrow, insert the batteries following the polarity as shown and close the lid.
Dispose of the used batteries as indicated in the paragraph "Disposal procedure".

⚠	When the batteries are completely depleted the appliance does not illuminate.
---	---

DATE/TIME SETTING

- 1) When the appliance is off, hold the SET button down until the display flashes the year. Use the "M" button to adjust the date and confirm with the "SET" button.
- 2) In succession, set the month, day, hours and minutes. Adjust the data with the "M" button and confirm with the "SET" button. The display will show the time in the 24 hour format.

When you replace the batteries, it is necessary to once again set the date/time settings

WRAPPING THE CUFF

- 1) Remove the cuff from its compartment.
- 2) Wrap the cuff on the left arm as shown in the figures "winding the cuff" and fasten the cuff with the tear-off closure.
The edge of the cuff should be about 2-3 cm above the elbow and the palm of your hand should be facing upward, resting on a surface.
Position the air hose at the centre of the cuff. The cuff should fit snugly around the arm but should not be pulled too tightly; leave enough room to insert a finger between the cuff and arm.

If the cuff is too tight or too loose, the readings will be inaccurate. Never roll sleeves up or the blood flow will be obstructed and the reading will not be accurate.

⚠	Make sure that the supplied cuff (22 - 42 cm) corresponds to the arm circumference.
---	---

PROPER READING METHOD

- For accurate blood pressure readings, follow these instructions:
- 1) Sit down, relax and keep still for at least 5 minutes before taking the reading.
 - 2) Remove sleeves and jewellery from the arm and wrist before applying the cuff.
 - 3) Avoid eating, smoking, drinking and exercising before taking the reading.
 - 4) Always use the same arm (preferably the left one) to take the reading.
Rest the arm on a hard surface so that the cuff is level with the heart.
The arm must be stretched out and relaxed.
When the reading is being taken, do not move any part of the body or the blood pressure monitor.
 - 5) Rest both feet on the floor without crossing legs or feet.
The reading can also be taken when lying down. It is sufficient to lie on your back, applying the cuff level with the heart, with the palm of your hand facing upwards (see the "correct reading position" figure).
- 6) Readings should always be taken at the same time to allow a comparison of blood pressure readings.**
- 7) Do not make reference to a single measurement. It is recommended to carry out at least 2 measurements, with intervals of at least 10/15 minutes from each other.**
It is necessary to let the arm rest for that amount of time since the congestion of blood may cause false readings. Do not perform more than 6 measurements per day.
- 8) Should you feel unpleasant sensations during a measurement, immediately turn off the device with the key "O/I".

TAKING A READING

- 1) Insert the plug of the air tube into the socket of the blood pressure monitor.
- 2) Press the "O/I" key. On the display, only for a few seconds, all the function symbols light up. The display shows "0". In the case in which air has remained in the cuff from the previous measurement, the symbol (🔻) flashes on the display for a few seconds.
- 3) The cuff inflates automatically and stops once it has reached the optimal level. Try

to stay relaxed, without speaking and without moving.
If the pre-set pressure (170 mmHg) is deemed insufficient or if there is a movement of the hand, the unit will carry out the inflation again (up to a maximum of 280 mmHg).

- 4) The cuff deflates automatically and on the display is shown the systolic pressure, diastolic pressure, pulsations, date and time of the measurement.
The symbol (♥) appears only in the case in which erratic heart beats (arrhythmias) have been detected. The segments that appear to the left of the display indicate the classification of the blood pressure values.
- 5) The device switches off automatically after approximately 1 minute of non-use or by pressing the "O/I" button.
It is possible to interrupt a measurement by pressing the "O/I" button.

⚠	Make sure that the batteries are fully charged: depleted batteries or those with little charge decrease the efficiency of the pump, which fails to give to the blood pressure monitor sufficient pressure for inflation in the pre-set time interval. For this reason the blood pressure monitor will indicate ERR. Therefore, replace the batteries.
---	--

MEMORY FUNCTION

This device is able to store up to 90 blood pressure measurements. After each measurement, all the detected values are automatically stored. To visualise past measurements, press the "M" button: the display shows the average of the measurements.
To scroll through the stored data, press the "M" or "SET" buttons and the measurements will be displayed, from the most recent to the oldest: the number 01 indicates the most recent data, the number 90 the oldest.
The device switches off automatically after approximately 1 minute of non-use or by pressing the "O/I" button. Having exceeded the 90 measurements, the oldest data are automatically cancelled. If there are no measurements stored, the display will show "n/a".

DELETION OF SAVED DATA

It is possible to erase all stored data: press the "M" button for about 5 seconds, the display will show "n/a" and it turns off automatically.

MAINTENANCE

- Store the device in its case in a well-ventilated and dry place without wrapping the hose too tightly and without resting any heavy objects on top of it.
- Clean the blood pressure monitor using a soft cloth, either dry or moistened with water and a liquid disinfectant.
- ⊗ Never use chemicals or abrasives.
- ⊗ NEVER wash the cuff in the washing machine and do not rub it vigorously, but rather delicately wipe the surface with a soft cloth, either dry or moistened with ethyl alcohol (75-90%) and leave it to dry naturally.
- Make sure liquids are never poured into the air hose.
- ⊗ Do not press the "O/I" key if the cuff is not wrapped around the arm.
- ⊗ Do NOT disassemble the device.
- It is recommended to check the performance of the device every 2 years. Contact Laica service assistance (activities excluded from the warranty).
- Repairs of the blood pressure monitor should not be made by the user.

TROUBLESHOOTING

Problem	Possible cause	Solution
After pressing the "O/I" key the reading does not start.	The batteries have not been inserted properly.	Check the correct positioning of the batteries.
	The batteries are worn.	Replace them.
	Strong electromagnetic interference.	Remove the batteries for 5 minutes and take the reading again.
The display shows the dead battery symbol  .	The batteries are worn. If normal zinc-carbon batteries are used, they will need to be replaced more often	Replace them. Use alkaline batteries to increase operating autonomy.
The readings are very low or very high.	The cuff is not correctly positioned.	Re-read the "Wrapping the cuff" paragraph.
	Incorrect posture while taking the reading.	Re-read the "Proper reading method" paragraph.

Problem	Possible cause	Solution
The readings are very low or very high.	During the reading, the person made a movement or spoke, or the reading was taken at a time when the person was particularly anxious or nervous.	Re-read the "Proper reading method" paragraph.
The heart rate readings are too low or too high.	A movement was made during the reading.	Re-read the "Proper reading method" paragraph.
	The reading was taken after exercise.	
The «♥» symbol appears on the display.	An irregular heartbeat was detected (arrhythmia).	Repeat the reading, if the symbol appears again, seek medical advice.
The display shows "E1".	The monitor is not properly inflating the cuff.	Check that the cuff is intact and that it is not leaking air. If damaged, buy a new cuff (part number ABM013).
The display shows "E3"	The cuff pressure has exceeded 299 mmHg.	Wait 5 minutes and repeat the measurement. In the case in which the display continues to display this error, please contact customer support.
The display shows "E2E4".	There is movement during the measurement.	Reread the paragraph "Correct method of measurement".

DISPOSAL PROCEDURE

 The symbol on the bottom of the device indicates the separated collection of electric and electronic equipment (**Dir. 2012/19/Eu-WEEE**).
 At the end of life of the device, do not dispose it as mixed solid municipal waste, but dispose it referring to a specific collection centre located in your

area or returning it to the distributor, when buying a new device of the same type to be used with the same functions. If the appliance to be disposed of is less than 25 cm, it can be returned to a retail location that is over 400 m² without having to purchase a new, similar device. This procedure of separated collection of electric and electronic devices is carried out forecasting a European environmental policy aiming at safeguarding, protecting and improving environment quality, as well as avoiding potential effects on human health due to the presence of hazardous substances in such equipment or to an improper use of the same or of parts of the same. Caution! The wrong disposal of electric and electronic equipment may involve sanctions.

To correctly dispose of batteries (**Dir. 2013/56/Eu**) **do not throw them into household waste** but rather dispose of them as special waste at specialised recycling centres. For more information about disposal of used batteries, contact the store where you purchased the equipment containing the batteries, the municipality or the local waste disposal service.

WARRANTY

This appliance is guaranteed for 2 years from the time of delivery of the goods, or another longer term envisaged by the national legislation of the consumer's residence. This provision complies with Italian and European legislation. The Laica products are designed for home use and must not be used in public venues. The warranty only covers manufacturing defects and does not apply if the damage is caused by an accidental event, incorrect use, negligence or misuse of the product. Use only the accessories supplied; the use of different accessories may result in invalidity of the warranty. Do not open the unit for any reason; in the case of opening or tampering, the warranty is definitively voided. This warranty does not apply to parts subject to wear or to the batteries when supplied. After 2 years from delivery, or another longer term envisaged by the national legislation of the consumer's residence, the warranty expires; in this case, the technical assistance interventions will be carried out against a fee. Information on technical assistance, whether under warranty or for a fee, can be requested by contacting info@laica.com. No payment will be due for repairs or replacements of products that fall within the terms of the warranty. In the event of faults, contact the retailer. Do NOT send the appliance directly to LAICA. All the operations under warranty (including those of replacement of the product or part thereof) will not prolong the duration of the original

period of warranty of the product replaced. The manufacturer declines any liability for any damage that may, directly or indirectly, be caused to persons, property or animals as a result of the non-observance of all the requirements established in the relevant instructions manual and concerning, especially, warnings relating to installation, use and maintenance of the appliance. Laica, in its constant commitment to improving its products, is entitled to changing without any notice, in whole or in part, its products in relation to production requirements, without this entailing any liability for Laica towards its dealers. For further information: www.laica.it.

STANDARDS

The product corresponds to the below standards: IEC 60601-1 (Medical electrical equipment - Part 1:General requirements for basic safety and essential performance), IEC 60601-1-2 (Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests), IEC 80601-2-30 (Medical electrical equipment – Part 2-30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers) .

 Manufactured by: Shenzhen Pango Medical Electronics Co., Ltd. No. 25 1st Industry Zone, Fenghuang Rd, Xikeng village, Henggang Town, Longgang District, Shenzhen, Guangdong China

EC	REP	Lotus NL B.V., Koningin Julianaplein 10, 1e Verd, 2595AA, The Hague, Netherlands
----	-----	--

Distributed by: Laica S.p.A.
Viale del lavoro, 10 - 36048 Barbarano Mossano (VI) – Italy
Phone: +39.0444.795314 - www.laica.com
Made in China

Guidance and manufacture's declaration electromagnetic emissions		
The DEVICE is intended for use in the electromagnetic environment specified below. The customer or the user of the DEVICE should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The DEVICE is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Conforms Class B	The DEVICE is suitable for use in all establishments domestic included and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Immunity aspects			
The DEVICE is intended for use in the electromagnetic environment specified below. The customer or the user of the DEVICE should assure that it is used in such an environment.			
Immunity test	Test level EN 60601-1-2	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) EN 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Burst/Fast Transient EN 61000-4-4	±2 kV for power supply lines and patient coupled lines	±2 kV for power supply lines and patient coupled lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge EN 61000-4-5	±1 kV line(s) and neutral	±1 kV line(s) and neutral	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines EN 61000-4-11	< 5% UT (>95% dip in UT) for 0,5 cycles 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles < 5% UT (>95% dip in UT) for 5 seconds	< 5% UT (>95% dip in UT) for 0,5 cycles 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles < 5% UT (>95% dip in UT) for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions. It is recommended that the DEVICE be powered from an uninterruptible power supply or a battery.
Power frequency magnetic field EN 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be that of a typical commercial or hospital environment
Note: UT is the a.c. mains voltage prior to application of the test level.			

Immunity aspects at r.f.			
The DEVICE is intended for use in the electromagnetic environment specified below. The customer or the user of the DEVICE should assure that it is used in such an environment.			
Immunity test	Test Level EN 60601-1-2	Compliance level	Electromagnetic environment – guidance
RF conducted EN 61000-4-6	3 Vrms from 150kHz to 80MHz	3 Vrms	Portable and mobile RF communications equipment shouldbe used no closer to any part of the device, including cables, than the recommended separation distance calculated from that equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1,2\sqrt{P}$ $d = 1,2\sqrt{P}$ 80 MHz to 800 MHz $d = 2,3\sqrt{P}$ 800 MHz to 2,5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).
RF radiated EN 61000-4-3	3 V/m from 80MHz to 2,5GHz	3 V/m	
Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, c should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: ((⦿))			

- NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.
- NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.
- Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the DEVICE is used exceeds the applicable RF compliance level above, the DEVICE should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the DEVICE.
 - Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the DEVICE			
The DEVICE is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the DEVICE can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the DEVICE as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1,2\sqrt{P}$	80 MHz to 800 MHz $d = 1,2\sqrt{P}$	800 MHz to 2,5 GHz $d = 2,3\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be determined using the equation to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer			

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NORMAS

El producto cumple con las siguientes normas: IEC 60601-1 (Medical electrical equipment - Part 1:General requirements for basic safety and essential performance), IEC 60601-1-2 (Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests), IEC 80601-2-30 (Medical electrical equipment – Part 2-30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers) .

 Producido por: Shenzhen Pango Medical Electronics Co., Ltd. No. 25 1st Industry Zone, Fenghuang Rd, Xikeng village, Henggang Town, Longgang District, Shenzhen, Guangdong China

  Lotus NL B.V., Koningin Julianaplein 10, 1e Verd, 2595AA, The Hague, Netherlands

Distribuido por: Laica S.p.A.
Viale del lavoro, 10 – 36048 Barbarano Mossano (VI) – Italy
Phone: +39.0444.795314 - www.laica.com
Made in China

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NORMAS

O producto atende às seguintes normas: IEC 60601-1 (Medical electrical equipment - Part 1:General requirements for basic safety and essential performance), IEC 60601-1-2 (Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests), IEC 80601-2-30 (Medical electrical equipment – Part 2-30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers) .

Produzido por: Shenzhen Pango Medical Electronics Co., Ltd. No. 25 1st Industry Zone, Fenghuang Rd, Xikeng village, Henggang Town, Longgang District, Shenzhen, Guangdong China

Lotus NL B.V., Koningin Julianaplein 10, 1e Verd, 2595AA, The Hague, Netherlands

Distribuído por: Laica S.p.A.
Viale del lavoro, 10 - 36048 Barbarano Mossano (VI) – Italy
Phone: +39.0444.795314 - www.laica.com
Made in China

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STANDARDS

Das Produkt erfüllt die folgenden Standards: IEC 60601-1 (Medical electrical equipment - Part 1:General requirements for basic safety and essential performance), IEC 60601-1-2 (Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests), IEC 80601-2-30 (Medical electrical equipment – Part 2-30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers).

 Hergestellt von: Shenzhen Pango Medical Electronics Co., Ltd. No. 25 1st Industry Zone, Fenghuang Rd, Xikeng village, Henggang Town, Longgang District, Shenzhen, Guangdong China

EC	REP
----	-----

Lotus NL B.V., Koningin Julianaplein 10, 1e Verd, 2595AA, The Hague, Netherlands

Vertrieben von: Laica S.p.A.
Viale del lavoro, 10 - 36048 Barbarano Mossano (VI) – Italy
Phone: +39.0444.795314 - www.laica.com
Made in China

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NORMES

Le produit répond aux normes suivantes: IEC 60601-1 (Medical electrical equipment - Part 1:General requirements for basic safety and essential performance), IEC 60601-1-2 (Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests), IEC 80601-2-30 (Medical electrical equipment – Part 2-30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers) .

 Produit par: Shenzhen Pango Medical Electronics Co., Ltd. No. 25 1st Industry Zone, Fenghuang Rd, Xikeng village, Henggang Town, Longgang District, Shenzhen, Guangdong China

EC	REP
----	-----

Lotus NL B.V., Koningin Julianaplein 10, 1e Verd, 2595AA, The Hague, Netherlands

Distribué par: Laica S.p.A.
Viale del lavoro, 10 - 36048 Barbarano Mossano (VI) – Italy
Phone: +39.0444.795314 - www.laica.com
Made in China

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Distribuirati: Laica S.p.A.