HC91/c - 03/2020

Cert. PG-IRT1603

LAICA S.p.A.
Viale del Lavoro, 10
36048 Barbarano Mossano (VI) - Italy
Phone +39 0444.795314
info@laica.com
Made in China

www.laica.it

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TH1004

TERMOMETRO A INFRAROSSI FRONTALE E AURICOLARE

EN	Instructions and warranty	page 12	RO	Instrucțiuni și garanție	pagina 60
	INFRARED EAR AND FOREHEAD THERMOMETER			TERMOMETRU CU INFRAROȘU DE URECHE ȘI FRUNTE	
ES	Instrucciones y garantía	página 20	CS	Návod k použití a záruční podmínky	strana 68
	TERMÓMETRO DE INFRARROJOS AURICULAR Y FRONTAL			INFRAČERVENÝ UŠNÍ A ČELOVÝ TEPLOMĚR	
PT	Instruções e garantia	página 28	SK	Návod na použitie a záruka	strana 76
	TERMÓMETRO INFRAVERMELHO AURICULAR E FRONTAL			INFRAČERVENÝ TEPLOMER NA MERANIE TEPLOTY V UCHU A NA ČELE	
DE	Anleitungen und garantie	Seite 36	HU	Utasítások és garancia	oldal 84
	INFRAROT - OHR- UND STIRNTHERMOMETER			INFRAVÖRÖS HOMLOK- ÉS NO CONTACT LÁZMÉRŐ	
FR	Instructions et garantie	page 44	SL	Navodila in garancija	Stran 93
	THERMOMÈTRE À INFRAROUGE AURICULAIRE ET FRONTALE			TERMOMETR NA INFRARDEČE ŽARKE ZA MERJENJE V UŠESU IN NA ČELU	
EL	Οδηγίες και εγγύηση	σελίδα 52	HR	Upute i jamstvo	Stranica 101
	ΘΕΡΜΟΜΕΤΡΟ ΥΠΕΡΥΘΡΩΝ ΑΥΤΙΟΥ ΚΑΙ ΜΕΤΩΠΟΥ			INFRAČRVENI TERMOMETAR ČEONI I UŠNI	

Fig. 1


Fig. 2


Fig. 3



Fig. 4

Fig. 5

LAICA[®] 2

ANNI DI **GARANZIA**
YEARS **GUARANTEE**
ANOS DE **ARANTIA**
ANOS DE **ARANTIA**
ANS DE **GARANTIE**
JAHRE **GARANTIE**
XPONIA **ΕΓΓΥΗΣΗ**
ANI DE **GARANȚIE**
ROK **ZÁRUKA**
ROK **ZÁRUKA**
ÉV **GARANCIA**
LET **GARANCIJE**
GODINA **GARANCIJE**

LAICA S.p.A.

Viale del Lavoro, 10
36048 Barbarano Mossano (VI) - Italy
Phone +39 0444.795314
info@laica.com
Made in China

www.laica.it

Data - Date

--	--	--	--	--	--

Timbro rivenditore
Retailer's stamp
Sello del revendedor
Carimbo do revendedor
Stempel des Händlers
Cachet du revendeur
Σφραγίδα καταστήματος

Ștampilă distribuitor
Razítko prodejce
Pečiatka predajcu
Bolti pecsét
Pečat prodajalca
Pečat prodavaoca

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INFRARED EAR AND FOREHEAD THERMOMETER
INSTRUCTIONS AND WARRANTY

Dear Customer, Laica would like to thank you for choosing this product, designed according to criteria of reliability and quality for complete satisfaction.

IMPORTANT
READ CAREFULLY BEFORE USING - KEEP IN A SAFE PLACE FOR FUTURE REFERENCE

The instruction manual must be considered as part of the product and must be kept for its entire life cycle. If the appliance is transferred to another owner, its documentation must also be transferred in its entirety. To ensure the safe and correct use of the product, the user must carefully read the instructions and warnings contained in the manual insofar as they provide important information regarding safety, user and maintenance instructions. In case of loss of the instruction manual or to receive further information or clarifications contact the company using the details below: LAICA S.p.A. - Viale del Lavoro, 10 - 36048 Barbarano Mossano (VI) - Italy - Phone +39 0444.795314 - info@laica.com - www.laica.it

Intended use and scope: The LAICA infrared thermometer detects the heat generated by the forehead and the ear canal to indicate the patient's body temperature. The innovative remote mode for front detection makes it ideal for measurements on children, even during sleep. It is very easy to use and its purpose is to detect body temperature and the temperature of liquids. It boasts a large LED display in four different colours according to the temperature detected. N.B.: the body temperature shown on the display is equivalent to the oral temperature.

KEY TO SYMBOLS

Warning

Attention! Read the user instructions carefully

Prohibition

Symbol of "type BF applied parts"

Manufacturer

Serial number

In compliance with European legislation on medical devices

Production batch number

IP22: Degree of protection of enclosures for electrical equipment, where the first digit indicates the degree of protection against the penetration of solid foreign bodies (from 0 to 6) and the second digit the degree of protection against the penetration of liquids (from 0 to 8).

SAFETY WARNINGS

- Before using the product, check that the appliance is intact without visible damage. If in doubt, do not use the product and contact your dealer.
- Keep the plastic bag away from children: danger of suffocation.
- This device must be used only for its intended purpose and according to the user instructions. All other types of use shall be considered inappropriate and therefore hazardous. The manufacturer shall not be deemed liable for any damage caused by inappropriate or incorrect use.
- People with physical, sensory or mental disabilities or unskilled people, may use and perform maintenance on the product only under adult supervision. Children must not play with the device.
- This appliance can be used autonomously by children of 12 years of age or older (able to read and understand the instructions for use). Below this age, they could damage their eardrum during ear measurement.
- Treat the product with care, protect it from impacts, extreme variations in temperature, humidity, dust, direct sunlight and sources of heat.
- In the event of a failure and/or malfunction, turn off the device without tampering with it. Always contact the retailer for repairs.
- Ensure hands are dry when using the on and off key. NEVER submerge the product in water or any other liquids.
- Device for domestic use that must not be used for first aid or for continuous temperature monitoring.

ATTENTION! PRIOR TO USING THIS DEVICE

FEVER

Fever should be considered a symptom and not a disease. In general, it is a sign that the body is fighting an infection. In fact, it is a defence mechanism of the body characterised by an increase in body temperature above 37°C and increased heart and breathing rates. Normal temperature within the body is about 37.5°C and this is in fact the value measured with the rectal method. The temperature measured in the mouth is about 0.5°C less (37°C) while the underarm recording is 1°C less (36.5°C). From this we can deduce the importance of a correct recording of the temperature that should be carried out by valid means, and if possible not at a time when the body temperature may be increased for physiological reasons (for example after meals).

In a healthy person the temperature is affected by many factors:

- individual character of the person (individual metabolism);
- age (infants and toddlers' body temperature is higher and decreases with age). In children, the temperature varies with greater intensity, speed and frequency;
- clothing;
- outside air temperature;
- time of day (in the morning body temperature is lower and increases during the day);
- motor and mental activities carried out;
- measurement method;
- phase of the menstrual cycle.

It should be noted that the thermometer and the person must have both been in the room, in which the measurement will take place, at least 30 minutes before taking the temperature. Room temperature must be between 10°C and 40°C. The table below shows a list of "normal" average temperatures: **it is however advisable to become familiar with recognising your normal temperature when you are physically well in order to be able to interpret the different values recorded.**

METHOD OF MEASURING AVERAGE NORMAL TEMPERATURE

UNDERARM	36,5°C
ORAL	37°C
RECTAL	37,5°C
EAR	37°C

- 1) **Self-measurement means checking, not diagnosis or treatment. Unusual values must always be discussed with your doctor. Under no circumstances should you alter the dosages of any medication prescribed by your doctor.**
- 2) This thermometer is a sensitive instrument. Handle with care and do not expose it to mechanical shocks. Do not squeeze, bend, drop or take apart the thermometer.
- 3) Do not touch with your fingers and do not blow on the probe; do not disassemble the thermometer.
- 4) Do not use this thermometer to detect rectal, oral or underarm temperature. It must be used to detect the forehead temperature at a distance of 0-5 mm, or the ear temperature.
- 5) Measurement accuracy can be compromised by: pronounced sweating of the forehead, taking vasoconstricting medication and skin irritations. Holding the thermometer in your hand for too long can cause a false recording of body temperature.
- 6) Do not use the device in the presence of ear infections, such as external otitis and tympanitis.
- 7) Clean the forehead where the probe will be pointed before each recording.
- 8) **Measure the temperature away from direct sunlight, wind and air conditioning.**
- 9) Do not measure temperature for thirty minutes after physical activity or after a bath.
- 10) The appliance may provide inaccurate readings when used in temperature or humidity conditions outside the limits indicated in the "Technical Specifications" paragraph.
- 11) Do not use near strong magnetic fields and keep away from radio systems or mobile phones (for more information on interferences see paragraph "Guidelines and declaration of the manufacturer").

SAFE USE OF BATTERIES

- Remove the batteries if the device is not used for long periods of time and keep them in a well-ventilated and dry place, at room temperature.
- ⚡ Do NOT use rechargeable batteries.
- ⚡ NEVER expose the batteries to sources of heat or direct sunlight. Failure to observe these instructions may result in damage to and/or the explosion of the batteries.
- ⚡ Do NOT throw the batteries into an open flame.
- Batteries must always be removed or replaced by adults only.
- Keep batteries out of reach of children: the ingestion of batteries can cause serious harm and death. If a battery is swallowed, seek urgent medical care.
- The acid in batteries is corrosive. Avoid contact with the skin, eyes or garments.

FOREHEAD 37°C

Each type of thermometer lends itself to the measurement of body temperature at a specific point of the body: the forehead thermometer only for forehead measurement; the ear thermometer only for measurement in the ear; the digital thermometer for traditional rectal, axillary and oral measurement.

Depending on the point where the measurement is made, a different temperature value is obtained. The variation can be between 0.2-1°C. You can not compare the temperatures measured with different measuring methods. It is therefore necessary to consider in case of self-diagnosis or to report to your doctor what type of thermometer and at which point of the body the temperature has been measured.

PRODUCT DESCRIPTION (see fig.1)

- 1) Measurement probe
- 2) Forehead measurement cap
- 3) Key "F"
- 4) "mem" key
- 5) Key "O"
- 6) Battery compartment
- 7) Batteries
- 8) LED display
 - a. Memory number
 - b. Acoustic signal indicator
 - c. Battery charge indicator
 - d. Indication of the measurement result
 - e. Unit of measure
 - f. Indication of liquid temperature measurement
 - g. Indication of forehead temperature measurement
 - h. Indication of ear measurement

TECHNICAL CHARACTERISTICS

- Product name: infrared thermometer

- Commercial name: TH1004
- Body temperature measurement range: from 34°C to 43°C (93.2°F - 109.4°F)
- Temperature and surface temperature measurement range: from 0°C to 93.2°C (32°F - 199.7°F)
- Display resolution: 0.1°C / 0.1°F
- LED display in 3 colours (for body temperature only):
 - 34-37.1°C / 93.2-98.8°F: green, normal temperature;
 - 37.2-38.1°C / 98.9-100.6°F: yellow, fever;
 - 38.2-43°C / 100.7-109.4°F: red, high fever
- Remote measurement: 0 – 5 mm
- Batteries: 2x1.5V alkaline AAA
- Switch off: automatic, after 30 seconds of non-use
- Battery life: approximately 1000 readings (turning the thermometer on and off after each measurement)
- Operating conditions: from 10°C to 40°C (50°F - 104°F), relative humidity 15%-93%, atmospheric pressure: 70 – 106 kPa
- Conditions of transportation and storage: from -25°C to 55°C (-13°F - 131°F), relative humidity ≤93%, atmospheric pressure: 50 – 106 kPa
- Tolerance:
 - from 35 to 42°C (95-107.6°F): +/- 0.2°C (0.4°F)
 - < 34.9°C (94.8°F) and > 42.1°C (107.7-°F): +/- 0.3°C (0.5°F)
- Forehead mode:
 - Reference body zone: forehead
 - Measurement zone: forehead
 - Clinical bias: -0.078
 - Acceptability limit: 0,243
 - Clinical repeatability: 0,069

INSTRUCTIONS FOR USE

INSERTION/REPLACEMENT OF THE BATTERY

This thermometer operates with 2 x 1.5V AAA alkaline batteries. When first used or when the display shows the battery symbol  or "LO", proceed with inserting or replacing

the batteries. Open the battery compartment by pressing the cover as shown in figure no.2. Insert the batteries bearing in mind the indicated polarity and close the cover. Dispose of used batteries as indicated in the "Disposal procedure" paragraph.

SETTINGS

• UNIT OF MEASUREMENT °C / °F

This thermometer can measure temperature in Celsius and Fahrenheit. When the thermometer is off keep the "mem" key pressed for approximately 6 seconds: the measurement unit will flash on the display (°C or °F). To change the unit of measurement, press the "mem" key. Then it is possible to take the measurement, or turn off the thermometer by pressing the key "O".

• SILENT DETECTION MODE

The thermometer is equipped with acoustic signal when the reading is complete. It is possible to set the silent detection mode making the thermometer ideal for measurements on children even during sleep.

Turn on the thermometer by pressing the "O" key, when the unit of measurement flashes on the display (°C/°F) press the "mem" key to remove the acoustic signal: "[]" disappears from the display. To activate the acoustic signal, press the "mem" key again when the unit of measurement (°C/°F) flashes on the display: the display shows "[]".

 When the silent detection mode is active the thermometer does not emit any acoustic signal after the measurements.

TEMPERATURE READING:

 The temperature shown on the display is the oral temperature.

The green LEDs indicate a normal temperature (between 34 and 37.1°C), the yellow colour indicates the presence of fever (temperature between 37.2 and 38.1°C), the red colour indicates high fever (with temperature between 38.2 and 43°C). The display shows the red "LO" symbols when the body temperature measured is below 34°C (93.2°F) and "HI" when it is above 43°C (109.4°F).

• FOREHEAD

Attention!

In order to ensure accurate recording of the temperature, clean the probe in its entirety with a cotton swab moistened with alcohol after each use. Make sure that the probe is completely dry before use.

AN INADEQUATE CLEANING OF THE PROBE CAN PROVIDE FALSE READINGS

- 1) Clean and dry the forehead where you will be positioning the probe.
- 2) Check that the forehead detection cap is inserted into the thermometer (see fig. 3).
- 3) Turn on the thermometer by pressing the "O" button, the thermometer will beep and the display will show the data of the last measurement taken, then emit a new acoustic signal and the display will flash: it is now possible to measure.
- 4) Point the probe at the centre of the forehead just above the space between the eyebrows at a distance of approximately 0-5 mm and press the "F" button. The device emits a beep and the display shows the measured value. Do not move the thermometer before the measurement performed beep. When the temperature exceeds 37.2°C the thermometer emits a series of short beeps.
- 5) It is possible to perform a new measurement when the LEDs change colour (white); the thermometer emits a beep and the unit of measurement (°C or °F) flashes.
- 6) The appliance switches itself off automatically after approximately 30 seconds or switch it off by pressing the "O" key

• EAR

- 1) Clean the ear where the probe is rested.
- 2) Remove the forehead measurement cap from the thermometer (see Fig. 4).
- 3) Switch on the thermometer by pressing the "O" key, the thermometer will emit a beep and the display will show the data of the last measurement and then emit a new beep and the display will flash: it is now possible to measure.
- 4) To obtain a precise measurement, gently pull the ear upwards to straighten the ear canal (see fig. 4) and adhere the probe to the latter aiming towards the latter membrane.
- 5) Press the "F" key: the appliance will beep and the display will show the measured value. Do not move the thermometer before the measurement performed beep. Remove the probe from the ear.
- 7) It is possible to perform a new measurement when the LEDs change colour (white); the thermometer emits a beep and the unit of measurement (°C or °F) flashes.
- 8) The appliance switches itself off automatically after approximately 30 seconds or switch it off by pressing the "O" key.

• **SURFACES AND LIQUIDS** 

This thermometer can record the temperature of warm surfaces or liquids (such as milk, water,...). It is possible to measure both with the forehead detection cap and without it.

- 1) When the thermometer is off, keep the “**On**” key pressed for approximately 6 seconds: the thermometer will sound and the display will show the last measurement, then emit a new beep and the display will flash: it is now possible to measure.
- 2) To measure the temperature, move the thermometer close to the relevant area at a distance of approximately 0-5 mm and press the “**0**” key. The device emits a beep and the display shows the measured value (white LED). Do not move the thermometer before the measurement performed beep.
- 3) It is possible to make a new measurement when the thermometer beeps and the unit of measurement (°C or °F) flashes.
- 4) The appliance switches itself off automatically after approximately 30 seconds or switch it off by pressing the “**On**” key.

 **Never immerse the probe in any kind of liquid. Always check the temperature of milk and of heated food before giving them to children and also the baby bath temperature.**

MEMORY FUNCTION

This appliance stores the last 9 measurements to assess any changes. Readings are automatically stored. Once 9 readings have been exceeded, the oldest data is automatically deleted.

To recall the stored measurements, when the appliance is off, press the “mem” key, each press of the key scrolls the data in memory: the number 1 refers to the most recent detection, 9 to the oldest detection. The appliance switches off automatically after approximately 30 seconds of non-use.

To delete the memorised data, when the device is off, press the “mem” key, then keep the “**On**” key pressed for approximately 5 seconds until the display shows “... ..”.

MAINTENANCE

After each use, clean the thermometer probe with a cotton swab soaked in alcohol to ensure precise measurements (see fig.5). The probe is the most delicate part of the thermometer: take particular care during the cleaning phases.

Allow the probe to dry completely for at least 30 minutes. The body of the thermometer should be cleaned with a soft, dry cloth. THE BODY OF THE THERMOMETER MUST NOT BE WASHED WITH WATER. Never use abrasive products and do not immerse the thermometer in water or other liquids. Keep the thermometer in a clean, dry place, away from direct sunlight.

PROBLEMI E SOLUZIONI		
Problem	Possible cause	Solution
The device does not turn on.	The batteries are worn.	Replace the batteries.
	The batteries have not been installed correctly.	Check the polarity of the batteries during insertion.
The display shows  .	The batteries are flat.	Replace the batteries.
The display shows “ Hi ” or “ Lo ”.	The measured temperature does not fall within the thermometer measurement range: • body measurement: from 34°C to 43°C, • surface and liquid measurement: from 0°C to 93.2°C.	
	The temperature of a surface or a liquid is being measured without changing the thermometer measurement setting.	For measurement of the temperature of surfaces or liquids it is necessary to select the detection  .
The display shows “ Lo ”.	During forehead detection, the measurement distance is greater than 5 mm.	Approach the thermometer to the forehead during measurement at a distance of between 0 and 5 mm.
	During ear measurement the probe is not correctly inserted into the ear (too far).	Position the thermometer correctly in the ear, read the paragraph "Ear temperature detection"  .

The display shows “ Lo ”.	The probe is dirty.	Clean the probe after each use as described in the "Maintenance" paragraph.
The display shows “ Er.H ” or “ Er.L ”.	The room temperature is not between 10°C and 40°C (50°F and 104°F).	Use the thermometer respecting the environmental conditions of operation (read the paragraph "Technical specifications").
The display shows “Err”	When measuring the temperature of surfaces or liquids, the thermometer is in an environment where the temperature varies rapidly up to 5°C.	The thermometer must be placed in the room where the measurement is taken at least 30 minutes before detection and placed in an area protected from fluctuations.
The measurements are low	During measurement the thermometer was positioned incorrectly.	Reread the paragraphs "Detection of the forehead temperature", "detection of the ear temperature".
	The probe or the ear canal is dirty.	Clean the probe after each use and clean the ear canal before measurement.

N.B. Please Note: When the device does not recover proper operation despite having performed these checks, contact the retailer.

DISPOSAL PROCEDURE (Dir. 2012/19/Eu-WEEE)

The symbol on the bottom of the device indicates the separated collection of electric and electronic equipment. At the end of life of the device, do not dispose it as mixed solid municipal waste, but dispose it referring to a specific collection centre located in your area or returning it to the distributor, when buying a new device of the same type to be used with the same functions.

If the appliance to be disposed of is less than 25 cm, it can be returned to a retail location that is over 400 m² without having to purchase a new, similar device. This procedure of separated collection of electric and

electronic devices is carried out forecasting a European environmental policy aiming at safeguarding, protecting and improving environment quality, as well as avoiding potential effects on human health due to the presence of hazardous substances in such equipment or to an improper use of the same or of parts of the same.

Caution! The wrong disposal of electric and electronic equipment may involve sanctions.To correctly dispose of batteries (Dir. 2013/56/Eu) do not throw them into household waste but rather dispose of them as special waste at specialised recycling centres. For more information about disposal of used batteries, contact the store where you purchased the equipment containing the batteries, the municipality or the local waste disposal service.

WARRANTY

This equipment is covered by a 2-year warranty starting from the purchase date, which has to be proven by means **of the stamp or signature of the dealer and by the fiscal payment receipt that shall be kept enclosed to this document.** This warranty period complies with the laws in force and applies only in case the consumer is an individual. Laica products are designed for domestic use only; therefore, their use for commercial purposes is forbidden. The warranty covers only production defects and does not apply, if the damage is caused by accidental events, misuse, negligence or improper use of the product.

Use the supplied accessories only; the use of other accessories invalidate the warranty cover. Do not open the equipment for any reason whatsoever; the opening or the tampering of the equipment, definitely invalidate the warranty. Warranty does not apply to parts subject to tear and wear along the use, as well as to batteries, when they are included. Warranty decays after 2 years from the purchase date. In this case, technical service interventions will be carried out against payment. Information on technical service interventions, under warranty terms or against payment, may be requested by writing at info@laica.com. Repair or replacement of products covered by the warranty terms are free of charge. In case of faults, refer to the dealer; DO NOT send anything directly to LAICA. All interventions covered by warranty (including product replacement or the replacement of any of its parts) shall not be meant to extend the duration of the original warranty period of the replaced product. The manufacturer cannot be made liable for any damages that may be caused, directly or indirectly, to people, things and pets resulting from the inobservance of the provisions contained in the suitable user manual, especially with reference to the instructions related to installation, use and maintenance of the equipment. Being constantly involved in improving its own products, Laica reserves the right to modify its product, totally or partially, with reference to any production need, without that this involves any responsibility whatsoever by Laica or its dealers.

Guidelines and declaration of the manufacturer

Manufacturer's declaration and guidelines - electromagnetic emissions		
The Device is intended for use in the electromagnetic environments specified below. The customer or user of the Device must make sure that the device is used in an appropriate environment.		
Emissions test	Compliance	Electromagnetic environment - guide
Radio frequency emissions CRISP 11	Group 1	The device only uses radio frequency energy for internal operation. Radio frequency emissions are therefore very low and do not cause any interference in electronic equipment located nearby.
Radio frequency emissions CRISP 11	Class B	The Device is used in home and it's powered by DC 3V
Harmonic emissions ICE 61000-3-2	Not applicable	
Voltage variations/flicker emissions ICE 61000-3-3	Not applicable	

Manufacturer's declaration and guidelines – electromagnetic Immunity			
The Device is intended for use in the electromagnetic environments specified below. The customer or user of the Device must make sure that the device is used in an appropriate environment.			
Immunity test	ICE 60601 test level	Level of compliance	Electromagnetic environment - guide
Electrostatic discharges (EDS) ICE 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors should be wood, concrete or ceramic. If the floors are covered with synthetic material, the relative humidity must be at least 30%.
Network frequency with magnetic fields (50 Hz/60 Hz) ICE 61000-4-8	30 A/m 50/60Hz	30 A/m 50/60Hz	Network frequency magnetic fields should have characteristic levels of a typical location in a commercial or hospital environment.
N.B. U _r is the mains voltage in a.c. before applying the test level.			

Manufacturer's declaration and guidelines – electromagnetic Immunity			
The Device is intended for use in the electromagnetic environments specified below. The customer or user of the Device must make sure that the device is used in an appropriate environment.			
Immunity test	ICE 60601 test level	Level of compliance	Electromagnetic environment - guide
Conducted radiofrequency IEC 61000-4-6	3 Vrms150 kHz to 80 MHz 6 Vrms 150 kHz to 80 MHz outside ISM bandsa	Not applicable	Portable and mobile radiofrequency communication devices must not be used at a distance "d" of any part of the Device, including cables, lower than the separation distance calculated from the equation relevant to the transmitter frequency. Recommended separation distance $d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$ $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$ 80 MHz to 800 MHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$ 800 MHz to 2,7 GHz
Conducted radiofrequency IEC 61000-4-3	10 V/m From 80 MHz to 2.7 GHz	10 V/m	Where P is the maximum transmitter output index in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed radiofrequency transmitters, as determined by an electromagnetic site survey ^a , must be lower than the compliance level in each frequency range ^b . Interference may occur near equipment marked with the following symbol:

NOTE 1 At 80MHz and 800MHz the maximum frequency range is applied.
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is influenced by the absorption and reflection of structures, objects and persons.
a. The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 3,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.
b. The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2,7 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 has been incorporated into the formulae used in calculating the recommended separation distance for transmitters in these frequency ranges.
c. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Model PG-IRT1603 Infrared Thermometer is used exceeds the applicable RF compliance level above, the Model PG-IRT1603 Infrared Thermometer should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Model PG-IRT1603 Infrared Thermometer.
d. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Separation distances recommended between portable and mobile RF communication devices and the Device.			
The device is intended for use in electromagnetic environments in which radiated radiofrequency disturbances are controlled. The customer or the user of the Device can contribute to the prevention of electromagnetic interference by maintaining a minimum distance between the portable and mobile radiofrequency communication devices (transmitters) and the Device based on the indications provided below, according to the maximum output power of communication devices.			
Maximum rated output power of the transmitter (W)	Separation distance depending on the frequency of the transmitter (m)		
	From 150 kHz to 80 MHz $d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$	From 80 kHz to 800 MHz $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$	From 800 MHz to 2,7 GHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
For transmitters whose maximum rated output power is not among those listed above, the recommended separation distance d in meters (m) can be determined by the equation relevant to the transmitter frequency, where P is the maximum rated output power of the transmitter, expressed in watts (W) according to the data provided by the transmitter manufacturer. NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is influenced by absorption and reflection of structures, objects and people.			

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